

Xeroform Gauze Dressing *Non-Adhering Fine Mesh Gauze*

Xeroform Gauze Dressing is composed of fine mesh gauze infused with 3% bismuth tribromophenate and petrolatum, offers dual-action wound care. Its antimicrobial properties prevent infections while maintaining a moist environment for rapid healing. Non-adherent and versatile, it's trusted by healthcare professionals for burns, abrasions, ulcers, and surgical wounds, ensuring patient comfort and promoting healing. Use beneath an absorbent appropriate secondary dressing that manages drainage and maintains a moist wound environment.



ENHANCED PATIENT COMFORT

Designed to allow the wound to breathe, minimizing skin and tissue maceration while maintaining a moist environment to encourage healing.



Not made with natural rubber latex.

Indications:

Xeroform Petrolatum Dressing is intended for use as a primary contact layer in dressing wounds such as lacerations, skin graft recipient sites, newly sutured wounds, abrasions, and minor or partial thickness burns. It may also be used as an initial layer in dressing surgical wounds with light exudate.

Features & Benefits:

- Impregnated with 3% bismuth tribromophenate and petrolatum
- Moist wound environment conducive to healing
- Non-adherent design reduces pain during dressing changes
- Clings and conforms to all body contours
- Sterile, single use device

Available Sizes:

Products	Size	Reorder #	Units
Xeroform Gauze Dressing	1in x 8in	PCHX3001	50/Box, 10 Boxes/Carton, 2 Cartons/Case
Xeroform Gauze Dressing	2in x 2in	PCHX3002	25/Box, 10 Boxes/Carton, 4 Cartons/Case
Xeroform Gauze Dressing	4in x 4in	PCHX3004	25/Box, 10 Boxes/Carton, 4 Cartons/Case
Xeroform Gauze Dressing	5in x 9in	PCHX3005	50/Box, 5 Boxes/Carton, 4 Cartons/Case



Compare to XEROFORM